



The **Aloha**® is a nasal-pillow style interface that comes with 3 different pillow sizes. This device is extremely lightweight and quiet, and offers an incredibly stable and comfortable seal. The **Aloha**® contains all of the components required to properly size the patient, eliminating the hassle of opening multiple products to get the right fit.



Product Features

- The ball-and-socket Elbow combined with Hi-Flex™ Tubing adapts to the most active of sleepers and accommodates various sleeping positions
- Exclusive Arced-Track™ technology allows the angle of the Pillow Reservoir to be adjusted for a more customized and natural fit
- Reinforced Headgear offers advantages over other headgears such as:
 - ◆ Over-the-Head strap can be worn in two different positions
 - ◆ Tubing Holder Ring allows for securing the tubing where you want it, such as over the head
 - ◆ Does not have uncomfortable “cheek flaps” or “non-breathing” silicone straps
 - ◆ Comes standard with Comfort Wraps™

Order Information

Part No.	Description	HCPCS Code	Code Description	HCPCS Frequency
ALO100	Aloha System – All Sizes Kit	A7034 + A7035	Nasal interface used with PAP device; Headgear used with PAP device	1 per 3 months + (see below)
ALO121	Aloha Replacement Pillows – SM	A7033	Pillows for nasal interface	2 pair per 1 month
ALO123	Aloha Replacement Pillows – MD			
ALO125	Aloha Replacement Pillows – LG			
ALO130	Aloha Replacement Headgear	A7035	Headgear used with PAP device	1 per 6 months
ALO131	Aloha Frame System – SM	A7034	Nasal interface used with PAP device	1 per 3 months
ALO133	Aloha Frame System – MD			
ALO135	Aloha Frame System – LG			

800-200-9842
www.innomedinc.com

InnoMed
TECHNOLOGIES

aloha™



The **Aloha**® Nasal-Pillow Interface is extremely lightweight and quiet, and offers an incredibly stable and comfortable seal. The Aloha offers the following key advantages over other interfaces:

Features	ALOHA	ResMed Swift LT	ResMed Swift FX	Respironics OptiLife	F&P Opus 360	F&P Pilairo
<i>No contact with the outside of the nose</i>	√	√	√	√	√	
<i>Pillow Reservoir can rotate to provide for nasal angle adjustments</i>	√	√				
<i>Has Ball-and-Socket Elbow</i>	√				√	√
<i>Has Hi-Flex™ Tubing</i>	√		√			√
<i>Headgear comes standard with Comfort Wraps™</i>	√		√			
<i>Weighs less than 65 grams</i>	√				√	√
<i>Headgear positioned away from the eyes</i>	√	√			√	√
<i>Tubing can be optionally secured over-the-head</i>	√	√				
<i>Features easy push in/pull out connection to Headgear</i>	√					
<i>Over-the-Head strap can be worn in two different positions</i>	√					
<i>Rotation axis is actually at the Pillows, not at the Reservoir</i>	√					

(Replacement Frame Systems, Pillows, and Headgears available)

Order : 800-200-9842
www.innomedinc.com

InnoMed
TECHNOLOGIES